
Z-CAD® One4All

| | | | |
|-------------|-----------|--------------------------------|----|
| EN | ENGLISH | Operating instructions:..... | 2 |
| DE | DEUTSCH | Gebrauchsanweisung:..... | 8 |
| FR | FRANÇAIS | Mode d'emploi:..... | 14 |
| IT | ITALIANO | Istruzioni per l'uso:..... | 20 |
| ES | ESPAÑOL | Instrucciones de uso:..... | 26 |
| PT | PORTUGUÊS | Manual de instruções:..... | 32 |
| RU | РУССКИЙ | Инструкция по применению:..... | 38 |
| ARAB | عربي | الاستخدام:..... | 44 |
| JP | 日本語バージョン | 取扱説明書:..... | 50 |
| CN | 中文 | 使用说明书圆片坯料:..... | 56 |

CE 0297

METOXIT AG
Emdwiesenstrasse 6
CH-8240 Thayngen
Switzerland
Info@metoxit.com
www.metoxit.com
SRN CH-MF-000019683

EC REP EC IMP

M+S Ceramics Kft.
Rozmaring utca 13.
HU-8000 Székesfehérvár
Hungary
Info@m-s-ceramics.com
SRN HU-AR-000006710 (EC-REP)
SRN HU-IM-000006711 (EC-IMP)

EN ENGLISH**Operating instructions****01. Product description**

Metoxit® **Z-CAD® One4All** is a zirconium dioxide ceramic material for the production of ceramic prosthetic restorations to restore the function and aesthetic of damaged or missing teeth.

02. Z-CAD® One4All products

Z-CAD® One4All
Z-CAD® One4All Multi
Z-CAD® One4All Multi-MX

03. Correct usage

Z-CAD® CAM blanks are designed for milling for the production of ceramic restorations to restore natural teeth or for placement on abutments.

Z-CAD® CAM blanks can be processed with CAD/CAM or copy milling systems. Z-CAD® zirconium oxide restorations offer excellent properties in terms of biocompatibility, stability and aesthetics. This makes them the first choice for producing crowns, inlays, onlays, veneers and bridges.

04. Indications

Z-CAD® One4All CAM blanks are pre-sintered, white or pre-colored and are intended for the fabrication of frameworks such as bars, bridges, crowns, inlays, onlays and veneers and meet class 5 according to ISO 6872:2019. Restorations milled from Z-CAD® one4All can also be stained with e.g. Z-CAD® Liquids for individualization or, after successful final sintering, be veneered e.g. with Prodigio Stain+Glaze Build Up Set.

Z-CAD® Deep Liquids are specially designed for the immersion process. Z-CAD® Surface Liquids are mainly used for individual staining and characterization.

Z-CAD® One4All Multi / Multi-MX CAM blanks are pre-sintered and intended for the production of fully anatomical or reduced crowns, veneers, inlays and onlays up to restorations with maximum spans.

With the characteristic strength of over 1000 MPa and a translucency value of 45%, the biocompatible high-performance zirconia corresponds to class 5 according to ISO 6872:2019.

Restorations milled from Z-CAD® one4All Multi / Multi-MX can also be stained with e.g. Z-CAD® Liquids for individualization or, after successful final sintering, be veneered with e.g. Prodigio Stain+Glaze Build Up Set.

Z-CAD® Deep Liquids are specially designed for the immersion process. Z-CAD® Surface Liquids are mainly used for individual staining and characterization.

For more information, see the relevant Z-CAD® Deep Liquid, Z-CAD® Surface Liquid and Prodigio Stain+Shade Build Up Instructions for Use at www.metoxit.com/company/download

05. Target patient group and intended users

Patients with a tooth defect or missing teeth, with erupted permanent dentition and adequate dental care, good health and no contraindications. The intended users are dentists/dental technicians who are thoroughly trained in dental restoration. Dental specialists must observe the local provisions and restrictions applicable to their profession.

Metoxit® Z-CAD® One4All are not intended for general use or as a commercially available product.

06. Contraindications

Allergies or hypersensitivities to chemical contents of the zirconium dioxide material (ZrO₂, Y₂O₃, HfO₂, Al₂O₃, SiO₂, Fe₂O₃, Na₂O, TiO₂, CaO, MgO, K₂O, Er₂O₃, Co₃O₄, MnO₂)

07. Warnings and precautionary measures

The prepared instructions are not sufficient as a sole means for processing and placing restorations. The restorations should only be processed and placed by dentists who are thoroughly trained in dental restoration. Processing and placement of restorations and associated components without appropriate training can lead to the failure and dislodgement of the restorations.

Take suitable precautionary measures to prevent aspiration and/or swallowing.

Wear gloves, mouth protection and safety glasses to avoid irritation. Observe the information in the current safety data sheet in the download section of Metoxit AG, www.metoxit.com/unternehmen/download

08. Caution

Z-CAD® One4All should not be used after the expiry date. See package labelling.

09. Residual risks and side effects

The clinical outcome of the dental treatment is influenced by several variables. The following possible residual risks and side effects can occur.

- Biting/chewing/prosthetic problems
- Bruising, damage to adjacent/opposite teeth
- Discomfort, gum injuries, local inflammations
- Local pain
- Recurrent cavities
- Local infection
- Local swelling
- Slight bleeding
- Longer recovery/healing time than expected
- Loss of prosthetic components
- Poor aesthetic result
- Recall in the dental practice
- Resorption of the upper/lower jawbone
- Risk of swallowing/inhaling small parts during the procedure

10. Chemical composition

| Elements | Unit | Z-CAD® One4All Z-CAD® One4All Multi Z-CAD® One4All Multi-MX |
|--------------------------------|------|---|
| ZrO ₂ | wt% | 90 – 95 |
| Y ₂ O ₃ | wt% | 6 – 8 |
| HfO ₂ | wt% | ≤ 5.0 |
| Al ₂ O ₃ | wt% | ≤ 0.5 |
| other oxides ¹ | wt% | ≤ 2.0 |

¹ other oxides; SiO₂, Fe₂O₃, Na₂O, TiO₂, CaO, MgO, K₂O, Er₂O₃, MnO₂, Co₃O₄

11. Physical properties

| Z-CAD® One4All Z-CAD® One4All Multi Z-CAD® One4All Multi-MX | |
|---|------------|
| Flexural strength (MPa) | ≥ 1000 |
| Density (g/cm ³) | ≥ 6.05 |
| CTE (10 ⁻⁶ K ⁻¹) | 10.5 ± 0.5 |

12. Indications for use

13. Indications for use Maximum number of pontics

| Z-CAD® One4All Z-CAD® One4All Multi Z-CAD® One4All Multi-MX indication | Construction* | total units | Connected pontics totally | Minimum wall thickness in mm | Connection cross-section in mm 2 |
|---|---------------|-------------|---------------------------|------------------------------|----------------------------------|
| Primary parts/double crowns | O | 1 | - | 0.5 | - |
| Single crown front tooth | O | 1 | - | 0.5 | - |
| Single crown posterior | O | 1 | - | 0.5 | - |
| anterior bridge | OXO | 3 | 1 | 0.5 | ≥ 9 |
| | OXXO | 4 | 2 | 0.7 | ≥ 12 |
| posterior bridge | OXO | 3 | 1 | 0.7 | ≥ 12 |
| | OXXO | 4 | 2 | 0.7 | ≥ 14 |

X indicated; digits show maximum number of pontics.

In Canada, bridge constructions are limited to 6 units with a maximum of 2 connected pontics.

Disclaimer:

The information relates to values in the densely sintered state.

If the recommended minimum parameters are undercut, there is no liability whatsoever for the manufactured structures.

The recommended minimum parameters are guide values that must be individually adjusted depending on the situation.

14. Indications for use Sintering instructions

The recommended sintering temperature is 1,500°C (Tab.1+2) **for all Z-CAD products.**

Speed sintering (Tab. 3) and a higher standard sintering temperature (Tab. 4) are possible.

Table 1 Individual objects

| Heating rate °C/min | T max °C | holding time min | cooling °C/min |
|---------------------|----------|------------------|------------------------------------|
| 10 | 1,500°C | 120 | Natural furnace cooling (10°C/min) |


Table 2 bridges

| Heating rate °C/min | T max °C | holding time min | cooling °C/min |
|---------------------|----------|------------------|-----------------------------------|
| 5 | 1,500°C | 120 | Natural furnace cooling (5°C/min) |


Table 3 Speed sintering

| Heating rate °C/min | T max °C | holding time min | cooling °C/min |
|---------------------|----------|------------------|----------------|
| 35 | 1,530°C | 30 | max. 45°C/min |

Individual objects and bridges up to 3 members;
≥4 links according to the dental technical situation


Table 4 Standard sintering temperature

| sinter range | Heating rate °C/min | T max °C | holding time min | cooling °C/min |
|--------------|---------------------|-------------------------|------------------|-------------------------|
| | tabel 1) + 2) | 1,500°C – 1,550°C | tabel 1) + 2) | Natural furnace cooling |

All Z-CAD® materials can be sintered between 1,500°C and 1,550°C. The requirements of ISO 6872:2019 for flexural strength are fully met. Depending on the furnace type selected, the translucency and color saturation can be individually influenced.


15. Indications for use General

- Frameworks for ceramic veneers should be constructed with cusp support and allow a uniform veneer layer thickness.
- Make sure that the wall thickness is sufficient and maximise the cross-sections of the connectors in a vertically-aligned oval shape.
- Observe the correct positioning and the appropriate milling strategy for the design and the blank in the CAM software.
- Insert the blank into the machine according to the operating instructions of the CAM system and enter the correct factor indicated on the collar of the blank.
- Start the machining process.
- Once the processing is finished, remove the blank with the manufactured framework from the CAM system and check it for damage.
- Use a tapered milling cutter to separate the framework.
- Remove any grinding dust with a clean brush and by carefully blowing it off with oil-free compressed air.
- After a colouring process has been carried out (follow the separate instructions), dry the framework before placing it in the sinter bowl.
- Place the frameworks for sintering in a sinter bowl filled with sinter balls and ensure that there is enough support for all restoration parts.
- Sinter large span bridges with sinter support structure in vertical position.
- Program the furnace according to its operating instructions and the Metoxit sintering table.
- Start the sintering process.
- Remove the framework from the oven only after it has cooled completely.
- Check the framework for defects, sufficient wall thickness and fit.
- If necessary, make minor adjustments with a water-cooled diamond tool.
- Rinse the framework under water and dry it.
- It is now ready for further processing.

16. Indications for use Colours

Restorations milled from Z-CAD® CAM blanks can be coloured with Metoxit Z-CAD® liquids. Z-CAD® Deep Liquids are specifically designed for the dipping method. Z-CAD® Surface Liquids are primarily used for individual colouring and characterisation. Further information is available in the relevant Z-CAD® Deep Liquid and Z-CAD® Surface Liquid operating instructions. www.metoxit.com/unternehmen/download

Illustration: Z-CAD® One4All Multi MX

This symbol on the edge of the Z-CAD® One4All Multi / Multi-MX indicates the colour gradient from light (top) to dark (bottom).

17. Transportation, Storage and Handling

Store the products in a dry place and at room temperature. Improper storage and transport can impair essential material properties and design features and lead to functional failure.

Before using for the first time, please check the packaging and contents for integrity. Check the declaration on the label and Z-CAD® One4All for equality. Never use damaged material. Always store Z-CAD® One4All in its original packaging, clean, dry and at room temperature, avoiding damage and contact with liquids.

18. Note to users and patients

Serious incidents related to the product must be reported to the manufacturer or its EC representative and to the competent authority of the Member State in which the user or patient is located.

19. Magnetic Resonance Imaging (MRI)

These products are manufactured from a material that is not affected by MRI and is MR safe. However, restorations produced from this material can be fitted with metal implants, abutments and screws that can be affected by MRI energy.

20. Further information

Additional information on Metoxit AG products can be found on the Metoxit website www.metoxit.com where relevant. Operating instructions can be downloaded from www.metoxit.com/unternehmen/download. Please note that milled dental restorations from Z-CAD® One4All CAM blanks are considered custom-made products.

21. Disposal

Disposal should be done in an environmentally friendly way in accordance with local regulations.

22. Information that must be provided to patients

Information on contraindications, warnings, preventative measures, side effects and complications with Metoxit products should be provided to patients. Please note that the dental restoration is milled. The patient must be informed on the MRI compatibility of the Metoxit products used.


23. Validity

On the publication of these operating instructions, all previous versions lose their validity.

24. Availability

Some Metoxit AG products are not available in all countries.

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26. Trademark

Z-CAD® is a registered name and trademark of Metoxit AG, Switzerland

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بين 1500 درجة مئوية و 1550 درجة مئوية © Z-CAD يمكن تلييد جميع مواد