
# Z-CAD<sup>®</sup> smile

<b>EN</b>	ENGLISH	Operating instructions: .....	2
<b>DE</b>	DEUTSCH	Gebrauchsanweisung: .....	8
<b>FR</b>	FRANÇAIS	Mode d'emploi: .....	14
<b>IT</b>	ITALIANO	Istruzioni per l'uso: .....	20
<b>ES</b>	ESPAÑOL	Instrucciones de uso: .....	26
<b>PT</b>	PORTUGUÊS	Manual de instruções: .....	32
<b>RU</b>	РУССКИЙ	Инструкция по применению: .....	38
<b>ARAB</b>	عربي	الاستخدام: .....	44
<b>JP</b>	日本語バージョン	取扱説明書: .....	50
<b>CN</b>	中文	使用说明书圆片坯料: .....	56

**CE** 0297

METOXIT AG  
Emdwiesenstrasse 6  
CH-8240 Thayngen  
Switzerland  
Info@metoxit.com  
www.metoxit.com  
SRN CH-MF-000019683

**EC REP** **EC IMP**

M+S Ceramics Kft.  
Rozmaring utca 13.  
HU-8000 Székesfehérvár  
Hungary  
Info@m-s-ceramics.com  
SRN HU-AR-000006710 (EC-REP)  
SRN HU-IM-000006711 (EC-IMP)

**EN ENGLISH****Operating instructions****1. Product description**

Metoxit® **Z-CAD® smile** is a zirconium dioxide ceramic material for the production of ceramic prosthetic restorations to restore the function and aesthetic of damaged or missing teeth.

**2. Z-CAD® smile products**

**Z-CAD® smile**  
**Z-CAD® smile Multi**  
**Z-CAD® smile Multi-MX**

**3. Correct usage**

Z-CAD® smile products are designed for milling for the production of ceramic restorations to restore natural teeth or for placement on abutments.

Z-CAD® smile can be processed with CAD/CAM or copy milling systems. Z-CAD® zirconium oxide restorations offer excellent properties in terms of biocompatibility, stability and aesthetics. This makes them the first choice for producing crowns, inlays, onlays, veneers and bridges.

**4. Indications**

**Z-CAD® smile** CAM blanks are pre-sintered, white or pre-colored and intended for the production of fully anatomical crowns, veneers, inlays and onlays up to three-unit bridges (side and front teeth). With the characteristic strength of over 610 MPa and a translucency value of 49%, the biocompatible high-performance zirconia corresponds to class 4 according to ISO 6872:2019.

Restorations milled from Z-CAD® smile can also be stained with e.g. Z-CAD® Liquids for individualization or, after successful final sintering, they can also be veneered with e.g. Prodigio Stain+Glaze Build Up Set.

**Z-CAD® smile Multi / Multi-MX** CAM blanks are pre-sintered and intended for the production of fully anatomical crowns, veneers, inlays and onlays up to three-unit bridges (side and front tooth area). With the characteristic strength of over 610 MPa and a translucency value of 49%, the biocompatible high-performance zirconia corresponds to class 4 according to ISO 6872:2019.

Restorations milled from Z-CAD® smile Multi / Multi-MX can also be colored with Z-CAD® Liquids for individualization, for example, or veneered with Prodigio Stain+Glaze Build Up Set, for example, after successful final sintering.

Z-CAD® Deep Liquids are specially designed for the immersion process. Z-CAD® Surface Liquids are mainly used for individual staining and characterization.

For more information, see the relevant Z-CAD® Deep Liquid, Z-CAD® Surface Liquid and Prodigio Stain+Shade Build Up Instructions for Use at [www.metoxit.com/company/download](http://www.metoxit.com/company/download)

**5. Target patient group and intended users**

Patients with a tooth defect or missing teeth, with erupted permanent dentition and adequate dental care, good health and no contraindications. The intended users are dentists/dental technicians who are thoroughly trained in dental restoration. Dental specialists must observe the local provisions and restrictions applicable to their profession.

Metoxit® Z-CAD® smile are not intended for general use or as a commercially available product.

**6. Contraindications**

Allergies or hypersensitivities to chemical contents of the zirconium dioxide material (ZrO<sub>2</sub>, Y<sub>2</sub>O<sub>3</sub>, HfO<sub>2</sub>, Al<sub>2</sub>O<sub>3</sub>, SiO<sub>2</sub>, Fe<sub>2</sub>O<sub>3</sub>, Na<sub>2</sub>O, TiO<sub>2</sub>, CaO, MgO, K<sub>2</sub>O, Er<sub>2</sub>O<sub>3</sub>, Co<sub>3</sub>O<sub>4</sub>, MnO<sub>2</sub>)

## 7. Warnings and precautionary measures

The prepared instructions are not sufficient as a sole means for processing and placing restorations. The restorations should only be processed and placed by dentists who are thoroughly trained in dental restoration. Processing and placement of restorations and associated components without appropriate training can lead to the failure and dislodgement of the restorations.

Take suitable precautionary measures to prevent aspiration and/or swallowing.

Wear gloves, mouth protection and safety glasses to avoid irritation. Observe the information in the current safety data sheet in the download section of Metoxit AG, [www.metoxit.com/unternehmen/download](http://www.metoxit.com/unternehmen/download)

## 8. Caution

Z-CAD® products should not be used after the expiry date, see packaging label.

## 9. Residual risks and side effects

The clinical outcome of the dental treatment is influenced by several variables. The following possible residual risks and side effects can occur.

- Biting/chewing/prosthetic problems
- Bruising, damage to adjacent/opposite teeth
- Discomfort, gum injuries, local inflammations
- Local pain
- Recurrent cavities
- Local infection
- Local swelling
- Slight bleeding
- Longer recovery/healing time than expected
- Loss of prosthetic components
- Poor aesthetic result
- Recall in the dental practice
- Resorption of the upper/lower jawbone
- Risk of swallowing/inhaling small parts during the procedure

## 10. Chemical composition

Elements	Unit	Z-CAD® smile Z-CAD® smile Multi Z-CAD® smile Multi-MX
ZrO <sub>2</sub>	wt%	88 – 93
Y <sub>2</sub> O <sub>3</sub>	wt%	8 – 10
HfO <sub>2</sub>	wt%	≤ 5.0
Al <sub>2</sub> O <sub>3</sub>	wt%	≤ 0.5
other oxides <sup>1</sup>	wt%	≤ 2.0

<sup>1</sup> other oxides; SiO<sub>2</sub>, Fe<sub>2</sub>O<sub>3</sub>, Na<sub>2</sub>O, TiO<sub>2</sub>, CaO, MgO, K<sub>2</sub>O, Er<sub>2</sub>O<sub>3</sub>, MnO<sub>2</sub>, Co<sub>3</sub>O<sub>4</sub>

## 11. Physical properties

Z-CAD® smile Z-CAD® smile Multi Z-CAD® smile Multi-MX	
Flexural strength (MPa)	≥ 600
Density (g/cm <sup>3</sup> )	≥ 6.05
CTE (10 <sup>-6</sup> K <sup>-1</sup> )	10.5 ± 0.5

**12. Indications for use**

**13. Indications for use Maximum number of pontics**

Z-CAD® smile Z-CAD® smile Multi Z-CAD® smile Multi-MX indication	Construction*	total units	Connected pontics totally	Minimum wall thickness in mm	Connection cross-section in mm <sup>2</sup>
Single crown front tooth	O	1	-	0.7	-
Single crown posterior	O	1	-	0.7	-
anterior bridge	OXO	3	1	0.8	≥ 14
	OXXO	-	-	-	-
posterior bridge	OXO	3	1	1	≥ 16
	OXXO	-	-	-	-

X indicated; digits show maximum number of pontics.

In Canada, bridge constructions are limited to 6 units with a maximum of 2 connected pontics.

Disclaimer:

The information relates to values in the densely sintered state.

If the recommended minimum parameters are undercut, there is no liability whatsoever for the manufactured structures.

The recommended minimum parameters are guide values that must be individually adjusted depending on the situation.


**14. Indications for use Sintering instructions**

The recommended sintering temperature is 1,500°C (Tab.1+2) **for all Z-CAD products.**

Speed sintering (Tab. 3) and a higher standard sintering temperature (Tab. 4) are possible.


**Table 1 Individual objects**

Heating rate °C/min	T max °C	holding time min	cooling °C/min
10	1,500°C	120	Natural furnace cooling (10°C/min)

**Table 2 bridges**


Heating rate °C/min	T max °C	holding time min	cooling °C/min
5	1,500°C	120	Natural furnace cooling (5°C/min)

**Table 3 Speed sintering**

Heating rate °C/min	T max °C	holding time min	cooling °C/min
35	1,530°C	30	max. 45°C/min

Individual objects and bridges up to 3 members;  
≥4 links according to the dental technical situation

**Table 4 Standard sintering temperature**

sinter range	Heating rate °C/min	T max °C	holding time min	cooling °C/min
	tabel 1) + 2)	1,500°C – 1,550°C	tabel 1) + 2)	Natural furnace cooling

All Z-CAD® materials can be sintered between 1,500°C and 1,550°C. The requirements of ISO 6872:2019 for flexural strength are fully met. Depending on the furnace type selected, the translucency and color saturation can be individually influenced.


### 15. Indications for use General

- Frameworks for ceramic veneers should be constructed with cusp support and allow a uniform veneer layer thickness.
- Make sure that the wall thickness is sufficient and maximise the cross-sections of the connectors in a vertically-aligned oval shape.
- Observe the correct positioning and the appropriate milling strategy for the design and the blank in the CAM software.
- Insert the blank into the machine according to the operating instructions of the CAM system and enter the correct factor indicated on the collar of the blank.
- Start the machining process.
- Once the processing is finished, remove the blank with the manufactured framework from the CAM system and check it for damage.
- Use a tapered milling cutter to separate the framework.
- Remove any grinding dust with a clean brush and by carefully blowing it off with oil-free compressed air.
- After a colouring process has been carried out (follow the separate instructions), dry the framework before placing it in the sinter bowl.
- Place the frameworks for sintering in a sinter bowl filled with sinter balls and ensure that there is enough support for all restoration parts.
- Sinter large span bridges with sinter support structure in vertical position.
- Program the furnace according to its operating instructions and the Metoxit sintering table.
- Start the sintering process.
- Remove the framework from the oven only after it has cooled completely.
- Check the framework for defects, sufficient wall thickness and fit.
- If necessary, make minor adjustments with a water-cooled diamond tool.
- Rinse the framework under water and dry it.
- It is now ready for further processing.

### 16. Indications for use Colours

Restorations milled from Z-CAD® smile can be coloured with Metoxit Z-CAD® liquids. Z-CAD® Deep Liquids are specifically designed for the dipping method. Z-CAD® Surface Liquids are primarily used for individual colouring and characterisation. Further information is available in the relevant Z-CAD® Deep Liquid and Z-CAD® Surface Liquid operating instructions.

**Illustration: Z-CAD® smile Multi-MX**

This symbol on the edge of the Z-CAD® smile Multi / Multi-MX indicates the colour gradient from light (top) to dark (bottom).

**17. Transportation, Storage and Handling.**

Store the products in a dry place and at room temperature. Improper storage and transport can impair essential material properties and design features and lead to functional failure.

Before using for the first time, please check the packaging and contents for integrity. Check the declaration on the label and Z-CAD® smile for equality. Never use damaged material. Always store Z-CAD® smile in its original packaging, clean, dry and at room temperature, avoiding damage and contact with liquids.

**18. Note to users and patients**

Serious incidents related to the product must be reported to the manufacturer or its EC representative and to the competent authority of the Member State in which the user or patient is located.

**19. Magnetic Resonance Imaging (MRI)**

These products are manufactured from a material that is not affected by MRI and is MR safe. However, restorations produced from this material can be fitted with metal implants, abutments and screws that can be affected by MRI energy.

**20. Further information**

Additional information on Metoxit AG products can be found on the Metoxit website [www.metoxit.com](http://www.metoxit.com) where relevant. Operating instructions can be downloaded from [www.metoxit.com/unternehmen/download](http://www.metoxit.com/unternehmen/download) Please note that milled dental restorations from Z-CAD® smile CAM blanks are considered custom-made products.

**21. Disposal**

Disposal should be done in an environmentally friendly way in accordance with local regulations.

## 22. Information that must be provided to patients

Information on contraindications, warnings, preventative measures, side effects and complications with Metoxit products should be provided to patients. Please note that the dental restoration is milled. The patient must be informed on the MRI compatibility of the Metoxit products used.

## 23. Validity

On the publication of these operating instructions, all previous versions lose their validity.

## 24. Availability

Some Metoxit AG products are not available in all countries.

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## 26. Trademark

**Z-CAD®** is a registered name and trademark of Metoxit AG, Switzerland

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بين 1500 درجة مئوية و 1550 درجة مئوية © Z-CAD يمكن تلبيد جميع مواد