

## Z-CAD<sup>®</sup> Liquid - The Intruder

The two for all options. Dipping or paint-on technique.



### Z-CAD<sup>®</sup> Deep liquid, Surface liquid

These ready to use, water soluble, infiltration liquids are specifically produced, to color the Metoxit AG white zirconia substructures. Regardless if dipping or paint-on methods are applied the coloring liquids emerge into the surface of unsintered substructures.

Z-CAD<sup>®</sup> Deep is primarily designed for the dipping technique of Z-CAD<sup>®</sup> HD and HTL blanks. Due to higher chroma values, the liquid gives the framework an ideal intensity for the veneering techniques and creates a core color from the depth. It is self-explanatory that the paint-on technique is also an option utilizing these liquids.

Z-CAD<sup>®</sup> Surface is primarily for the brush technique (paint-on) of the Z-CAD<sup>®</sup> One4All and *smile* blanks intended. The higher translucency value causes a lower chroma for monolithic restorations to set accents onto the surface. It is self-explanatory that the dipping technique is also an option utilizing these liquids.

## Z-CAD® Deep liquid Dentin

Art. No.	Product Deep liquid A1 - D4	Content
411950	Z-CAD Deep liquid A1	100ml
411951	Z-CAD Deep liquid A2	100ml
411952	Z-CAD Deep liquid A3	100ml
411953	Z-CAD Deep liquid A3.5	100ml
411954	Z-CAD Deep liquid A4	100ml
411955	Z-CAD Deep liquid B1	100ml
411956	Z-CAD Deep liquid B2	100ml
411957	Z-CAD Deep liquid B3	100ml
411958	Z-CAD Deep liquid B4	100ml
411959	Z-CAD Deep liquid C1	100ml
411960	Z-CAD Deep liquid C2	100ml
411961	Z-CAD Deep liquid C3	100ml
411962	Z-CAD Deep liquid C4	100ml
411963	Z-CAD Deep liquid D2	100ml
411964	Z-CAD Deep liquid D3	100ml
411965	Z-CAD Deep liquid D4	100ml

## Z-CAD® Deep liquid Special

Art. No.	Product Deep liquid Special	Content
411966	Z-CAD Deep liquid blue1	100ml
411967	Z-CAD Deep liquid blue2	100ml
411968	Z-CAD Deep liquid grey1	100ml
411969	Z-CAD Deep liquid grey2	100ml
411970	Z-CAD Deep liquid yellow	100ml
411971	Z-CAD Deep liquid orange	100ml
411972	Z-CAD Deep liquid brown	100ml
411973	Z-CAD Deep liquid pink	100ml



## Z-CAD® Deep liquid Set

Art. No.	Product Deep liquid Set	Content
411974	Z-CAD Deep liquid dentin set	16x 100ml
411975	Z-CAD Deep liquid effect set	8x 100ml

## Z-CAD® Surface liquid dentin

Art. No.	Product Surface liquid A1 - D4	Content
411976	Z-CAD Surface liquid A1	100 ml
411977	Z-CAD Surface liquid A2	100 ml
411978	Z-CAD Surface liquid A3	100 ml
411979	Z-CAD Surface liquid A3.5	100 ml
411980	Z-CAD Surface liquid A4	100 ml
411981	Z-CAD Surface liquid B1	100 ml
411982	Z-CAD Surface liquid B2	100 ml
411983	Z-CAD Surface liquid B3	100 ml
411984	Z-CAD Surface liquid B4	100 ml
411985	Z-CAD Surface liquid C1	100 ml
411986	Z-CAD Surface liquid C2	100 ml
411988	Z-CAD Surface liquid C3	100 ml
411989	Z-CAD Surface liquid C4	100 ml
411990	Z-CAD Surface liquid D2	100 ml
412002	Z-CAD Surface liquid D3	100 ml
411991	Z-CAD Surface liquid D4	100 ml

## Z-CAD® Surface liquid effect

Art. No.	Product Surface liquid effect	Content
411992	Z-CAD Surface liquid blue1	100 ml
411993	Z-CAD Surface liquid blue2	100 ml
411994	Z-CAD Surface liquid grey1	100 ml
411995	Z-CAD Surface liquid grey2	100 ml
411996	Z-CAD Surface liquid yellow	100 ml
411997	Z-CAD Surface liquid orange	100 ml
411998	Z-CAD Surface liquid brown	100 ml
411999	Z-CAD Surface liquid pink	100 ml



## Z-CAD® Surface liquid set

Art. No.	Product Surface liquid set	Content
412000	Z-CAD Surface liquid dentin set	16 x 100 ml
412001	Z-CAD Surface liquid effect set	8 x 100 ml

This perfectly tuned system of Z-CAD® Deep and Surface liquids offers easy handling and a high reproducibility for to be colored restorations. The Z-CAD® DEEP and Surface liquid Set is comprised of the 16 dentin colors, according to the \*Vita classical A1- D4 color system. We also offer a Z-CAD® Deep and Surface liquid "effect set" which is comprised of eight colors for individualization.

\*VITA classical A1 - D4® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG in Bad Säckingen.